

**ANALYTICAL PROCEDURE**  
**DRY PREMIX FOR WHITE-COLOUR COATING 37781 RBC**  
**(E000658-013)**

*The methods and reagents R referred to in this technique are described in the European Pharmacopoeia.*

Mixture of coating excipients to be dispersed in water before use.

Percentage composition (% *m/m*) is the following:

Glycerol	4.50
Hypromellose	74.80
Macrogol 6000	1.80
Magnesium stearate	4.50
Titanium dioxide (E 171)	14.40

## **1. CHARACTERS**

White coloured granules or powder.

## **2. IDENTIFICATION**

Dissolve the residue obtained during the "Total ash" test in a mixture of 8 ml of concentrated sulfuric acid *R* and 2 ml phosphoric acid *R*. Heat until the mixture boils. Centrifuge or filter. Pour 1 ml of the clear solution into 5 ml of water to which is added several drops of a solution of hydrogen peroxide *R*. A yellow to yellow-orange coloration develops.

## **3. TEST**

### **3.1. Suspension S**

Prepare a bubble-free homogeneous aqueous suspension containing 15 per cent *m/m* of the substance to be examined.

### **3.2. Viscosity (2.2.10)**

300 mPa.s to 900 mPa.s, determined on the suspension S with a rotating viscometer (Brookfield-type LVT) at 20 °C.

### **3.3. Dispersion**

Allow 665 g of the suspension S to pass through a tared sieve with a mesh size of 150 µm. Rinse several times with water then dry and weigh the sieve. The mass of the residue on the sieve does not exceed 0.5 g (0.5 per cent *m/m*).

### **3.4. Film formation**

Spread out part of the suspension used during the "[Dispersion](#)" test on a metallic or glass plate. A white film is formed after drying.

### **3.5. Loss on drying (2.2.32)**

Not more than 6.0 per cent, determined on 1.5 g by drying in an oven at 100 °C to 105 °C for 3 h.

### **3.6. Total ash (2.4.16)**

12.0 per cent to 17.0 per cent, determined on 1.5 g at a maximum temperature of 800 °C ± 25 °C.

### **3.7. Microbial contamination (If tested)**

Total aerobic microbial count (TAMC) (2.6.12) not more than  $10^3$  CFU/g and total combined yeasts/moulds count (TYMC) not more than  $10^2$  CFU/g. It complies with the test for *Escherichia coli* (2.6.13).

## **4. STORAGE**

Store in a well-closed container.